




# Understanding the Report

## Com |PI|i|t DX®

**On the first page of the report**  
on the results section

You will find all the genes  
that were analyzed and the genomic alterations,  
identified on the patient's tumor.


Results	
<b>Tumor Type: Melanoma</b>	
Genomic Alterations Identified	
<b>BRAF c.1799TA (p.V600E)</b> <b>KIT c.1900CT (p.R634W)</b> <b>P53 c.380CT (p.S127F), c.568 CT (p.P190S), c625AG (p.R209G)</b>	
Genomic Alterations Identified	
<b>ABL1, EZH2, JAK3, PTEN, AKT1, FBXW7, IDH2, PTPN11, ALK, FGFR1, KDR, RB1, APC, FGFR2, KIT, RET, ATM, FGFR3, KRAS, SMAD4, BRAF, FLT3, MET, SMARCB1, CDH1, GNA11, ML1, SMO, CDKN2A, GNAS, MPL, SRC, CSF1R, GNAQ, NOTCH1, STK11, CTNNB1, HNF1, NPM1, TP53, EGFR, HRAS, NRAS, VHL, ERBB2, IDH1, PDGFRA, ERBB4, JAK2, PIK3CA, ALK fusions, ROS1 fusions, RET fusions, NTRK1 fusions.</b>	



**On the second page of the report**  
on the interpretation section

You will find the genomic alterations  
with the respective on-and off-label approved therapies,  
and the respective Clinical trials, where available.

Genomic Alteration	Approved therapies in	Approved therapies in	Potential
<b>BRAF p.V600E</b>	ABRAFENIB, TRAMETINIB, VEMURAFENIB	Rigorafenib	Yes, see clinical trials section
<b>KIT p.R634W</b>	None	Imatinib, Regorafenib, Axitinib, Cabozantinib, Sorafenib, Nilotinib, Sunitinib	Yes, see clinical trials section
<b>P53 p.S127F, p.P190S, p.R209G</b>	None	None	None





# Understanding the Report

## Com | P I | i | t D X<sup>®</sup>

**Com|PI|i|t DX**  
Comprehensive Panel for Individualized Treatment

**GENEKOR**  
Committed to Biotechnological Innovation

Name: \_\_\_\_\_ Date Received: \_\_\_\_\_  
Sample's Code: \_\_\_\_\_ Date of Report: \_\_\_\_\_  
Report No: \_\_\_\_\_

**Genomic Alterations**

**BRAF Mutations**  
Somatic mutations in BRAF have been found in 87.5% of all melanoid melanomas (COSMIC; Davies et al. 2000; Hoels et al. 2012; Knaflhammer et al. 2012; Malabaroto et al. 2008). BRAF mutations are found in all melanoma subtypes but are the most common in melanomas derived from skin without chronic sun-induced damage (Curtin et al. 2005; Malabaroto et al. 2008). In this category of melanoma, BRAF mutations are found in 79% of samples (Curtin et al. 2005).  
The most prevalent BRAF mutations detected in melanoma are missense mutations that introduce an amino acid substitution of valine 600. Approximately 80-90% of V600 BRAF mutations are V600E (valine to glutamic acid; COSMIC; Lowy et al. 2012; Subraman et al. 2010). BRAF V600E mutations are associated with increased sensitivity to BRAF inhibitors (Chapman et al. 2011; Falchook et al. 2012a; Fisher et al. 2010; Fisher et al. 2012a; Hauschild et al. 2012; Sosman et al. 2012). Patients whose tumors harbored V600E and V600K mutations showed better responses to the BRAF inhibitor, vemurafenib, than to chemotherapy (docetaxel or paclitaxel; Fisher et al. 2012b). Patients with V600E or V600K-mutated tumors also showed better responses to trametinib than patients with BRAF wild type tumors (Falchook et al. 2012b).  
The FDA has approved use of the combination therapy dabrafenib and trametinib for melanoma patients with V600E or V600K mutations based on an interim analysis of a phase III trial showing improved response rates and response duration compared to dabrafenib alone (FDA 2014; Long et al. 2014). Improved progression-free survival and response rates compared to dabrafenib monotherapy were observed in an earlier phase I trial, as well (Fisher et al. 2012a). Patients with V600-mutated tumors and no prior BRAF inhibitors treated with combination vemurafenib and cobimetinib experienced improved response rate and progression-free survival compared to patients with V600-mutated tumors and recent progression on vemurafenib (Bos et al. 2014). This study included patients with V600E and V600K mutations (Bos et al. 2014).

**BRAF Mutations**  
GENEKOR Medical S.A. | 52 Spaton Ave. Gerakas, 15344 Athens | info@genekor.com www.genekor.com |  
Tel. (+30) 210 6032138 Fax. (+30) 210 6032148  
A company Certified and Accredited by EN ISO 9001:2008. Certificate Number: 04188204 and by ISO 13188. Certificate Number: 822

**Com|PI|i|t DX**  
Comprehensive Panel for Individualized Treatment

**GENEKOR**  
Committed to Biotechnological Innovation

Name: \_\_\_\_\_ Date Received: \_\_\_\_\_  
Sample's Code: \_\_\_\_\_ Date of Report: \_\_\_\_\_  
Report No: \_\_\_\_\_

(12/24) TP-53 mutations reported. Interestingly, more than half of these (N=41) were UV signature changes (Hecker and Iaco 2007) suggesting that TP53 is, in fact, a melanoma-promoting target of UV.

**Clinical Trials to Consider**

**IMPROVANI** The information indicated below reflects the information that is available at the time of this report. Please note that this information is constantly updated, so please visit the above mentioned website in order to get updated information on available clinical trials.

**BRAF-Associated Clinical Trials**

Trial Code	Phase	Trial Title
NCT01726138	II	COCC-1128: Open Label Phase I Trial of the BRAF Inhibitor Dabrafenib and the MEK Inhibitor Trametinib in Inoperable Stage II and Stage IV BRAF-Mutant Melanoma: Correlation of Response With the Kinome and Functional Mutations
NCT01939396	II	Circulating Melanoma Cells in Melanoma: Patients Treated with Selective BRAF Inhibitors
NCT01924603	II	Study Comparing Combination of Dabrafenib Plus Trametinib Versus Vemurafenib and Dabrafenib Monotherapy in BRAF-Mutant Melanoma
NCT02034026	II	Study of Neoadjuvant Use of Vemurafenib Plus Combination for BRAF-Mutant Melanoma with Tumor Lymph Node Metastases
NCT01792326	II	Dabrafenib, Trametinib, BRAF or BRAF/MEK Inhibitors with BRAF and BRAF V600E Mutated Melanoma with Brain Metastases: Biomarkers and Mechanisms
NCT01972347	II	Neoadjuvant Dabrafenib Plus Trametinib in AJCC Stage IIIc BRAF V600E Melanoma: Correlative Study
NCT01842527	II	Utility of Next-Gen BRAF Test of Melanoma
NCT01414132	II	Access Study of Trametinib for Subjects with Advanced Unresectable Stage IIIc or Stage IV Melanoma: Stage III BRAF V600E Mutation Positive Cutaneous Melanoma
NCT02009181	II	Dabrafenib and Trametinib in Treating Patients with Stage IV BRAF-Mutant Melanoma That Cannot be Removed by Surgery
NCT02009181	II	Study to Evaluate Treatment of Advanced-Stage Melanoma in Subjects with BRAF-Mutant Cutaneous Melanoma That Can Be Removed by This Study

GENEKOR Medical S.A. | 52 Spaton Ave. Gerakas, 15344 Athens | info@genekor.com www.genekor.com |  
Tel. (+30) 210 6032138 Fax. (+30) 210 6032148  
A company Certified and Accredited by EN ISO 9001:2008. Certificate Number: 04188204 and by ISO 13188. Certificate Number: 822

**Com|PI|i|t DX**  
Comprehensive Panel for Individualized Treatment

**GENEKOR**  
Committed to Biotechnological Innovation

Name: \_\_\_\_\_ Date Received: \_\_\_\_\_  
Sample's Code: \_\_\_\_\_ Date of Report: \_\_\_\_\_  
Report No: \_\_\_\_\_

**Genes Analyzed**

22 gene alterations													
ABL1	EGFR	JAK3	PTEN	ARX1	RSK2	SH2	PIK3R1	ALK	FGFR3	KIF5B	BR	APC	FGFR2
RET	MEF	ARF	FGFR4	KIF5B	SMAD4	BRG1	SMAD4	BRG1	SMAD4	BRG1	SMAD4	BRG1	SMAD4
PTEN	MEF	SMAD4	CDMT1	SMAD4	BRG1	SMAD4	BRG1	SMAD4	BRG1	SMAD4	BRG1	SMAD4	BRG1
CDKN2A	GNAS	MAPK	SHC	CSF1R	GNAS	GNAS	GNAS	GNAS	GNAS	GNAS	GNAS	GNAS	GNAS
ERBB1	CTNNA1	HR23A	HR23B	HR23C	HR23D	HR23E	HR23F	HR23G	HR23H	HR23I	HR23J	HR23K	HR23L
NRAS	SH2	ERBB2	IKK1	PDGFRA	ERBB4	JAK2							
PIK3CA													

**4 fusion transcripts**

ALK	ROS1	RET	NTRK1
-----	------	-----	-------

**Methodology**

The material for this analysis is the patient's Formalin-Fixed Paraffin Embedded (FFPE) tumor tissue.  
DNA was extracted from the sample under investigation using the QIAamp DNA FFPE Tissue kit. RNA was extracted using the High Pure RNA Micro Kit (Roche Life Science).  
Mutation hotspot regions of 50 genes were amplified using Ion AmpliSeq™ HotSpot Cancer Research Panel (Life Technologies). Copy number variations, SNPs, and indels were analyzed. Additionally, ALK, ROS1, RET and NTRK1 fusions expression were tested using Ion AmpliSeq™ RNA Fusion Lung Cancer Research Panel (Life Technologies) and FISH and CISH. Sequencing was carried out using the next Generation Sequencing platform Ion proton (Life Technologies). The sensitivity of the method is 5% of mutant allelic content.

GENEKOR Medical S.A. | 52 Spaton Ave. Gerakas, 15344 Athens | info@genekor.com www.genekor.com |  
Tel. (+30) 210 6032138 Fax. (+30) 210 6032148  
A company Certified and Accredited by EN ISO 9001:2008. Certificate Number: 04188204 and by ISO 13188. Certificate Number: 822

On the **Genomic Alterations Section** a more detailed interpretation about the Genomic Alterations identified is contained, including frequency ratios, biological impact and related targeted therapies.

On the **Clinical Trials Section** all the currently available clinical trials are reported according to the molecular profile of the patient's tumor.

Furthermore, details about the methodology of the assay are reported on the **Methodology Section**.