

## Specimen Preparation Instructions



**GeneKor**

Committed to Biotechnological Innovation

The key sample consideration for GeneKor is to submit a total mass of cells that is sufficient to extract the amount of DNA necessary for analysis. This checklist was designed to help us work together to identify the best sample for each patient. However, solid tumor samples vary greatly from patient to patient, and you may have unique challenges not covered here. Please contact our Client Service Representative with additional questions.

### SAMPLE TYPE

#### BLOOD SAMPLE

Blood should be collected in EDTA tubes. A minimum of two vials are necessary for testing. Blood collection can be performed at any qualified center of your choice. GeneKor does not perform blood collection at its premises or in any other location.



### BLOOD SAMPLE STORAGE

The blood sample can be stored at 4°C in your home refrigerator.

\* **DO NOT** freeze the sample.



**Do not freeze**

### SHIPPING INSTRUCTIONS

1. Place the GeneKor Requisition Form and all relevant patient insurance materials and reports into the GeneKor Specimen Kit.

2. Place the GeneKor Specimen Kit (including the GeneKor Requisition Form and relevant patient reports and insurance documents) into the provided shipping pack and seal the shipping pack.

3. Ship completed specimen kit (requisition form, insurance information and samples) to:

**GeneKor Medical S.A.  
52 Spaton Ave.  
153 44 Gerakas  
Athens, Greece**

Place all materials into GeneKor Specimen Kit and contact Client Service to schedule a pick up.

[www.genekor.com](http://www.genekor.com)

Please call **+30 210 603 2138** with questions or to request GeneKor Specimen Kits.

[info@genekor.com](mailto:info@genekor.com)