

Specimen Preparation Instructions



Genekor
Committed to Biotechnological Innovation

This checklist was designed to help us work together to identify the best sample for each patient. Please contact our Client Service Representative if you have additional questions.

SAMPLE TYPE

BLOOD SAMPLE

Blood should be collected in EDTA tubes. A minimum of two vials are necessary for testing. Blood collection can be performed at any qualified center of your choice. Genekor does not perform blood collection at its premises or in any other location.



BLOOD SAMPLE STORAGE

The blood sample can be stored at 4°C in your home refrigerator.

* **DO NOT** freeze the sample.



Do not freeze

SHIPPING INSTRUCTIONS

* During the summer period we recommend to have an ice pack in the Kit (The ice pack **SHOULD NOT** directly touch the sample).

1. Place the GeneKor Requisition Form and all relevant patient insurance materials and reports into the GeneKor Specimen Kit.

2. Place the GeneKor Specimen Kit (including the GeneKor Requisition Form and relevant patient reports and insurance documents) into the provided shipping pack and seal the shipping pack.

3. Ship completed specimen kit (requisition form, insurance information and samples) to:

Genekor Medical S.A.
52 Spaton Ave.
153 44 Gerakas
Athens, Greece

Place all materials into GeneKor Specimen Kit and contact Client Service to schedule a pick up.

* **WARNING, vials that are not labeled with the full name and/or code of the patient will NOT be accepted.**

www.genekor.com

Please call **+30 210 603 2138** with questions or to request GeneKor Specimen Kits.

info@genekor.com